



Specifications

Pump Type	Infusion
Max.No. Syringes	One
Syringe Size	.5µL to 60mL
Linear Force	22lbs.
Min Flow Rate	> .1µL/hr
Max Flow Rate	8.87mL/min
Accuracy	± 3%
Reproducibility	± 1%
Power	9Volt DC @ 1000mA
Dimensions	9.0" x 5.75" x 4.93"
Weight	3.9 lbs.
Steps per Revolution	400

IN-1000 Advance Infusion Pump

The **Advance Infusion Pump Series 1200** is a simple, accurate, affordable single syringe pump. Easily input the syringe diameter, volume, and time into the large number keypad that no other pump on the market has and let Advance do the rest.

No need to scroll through long lists of syringe sizes and infusion rates again with Advance!

IN-1010 Advance Infusion Table

The **Advance Infusion Table** is the latest design in tail infusions for rats. This unique setup has allowed for tail infusions lasting up to three hours.

Features

- 4 folding legs for optional angled display and compact storage
- Easy customization using the silicone tubing and perforated table top
- Cage washable

Dimensions

15 3/4" x 15 3/4" square table top (40cm x 40cm)
 Height with legs open, 6 1/2" (16.5cm)

