

RS-9900	Micro Sterilizing 8" x 8 ^{1/2} " x 1 ^{1/4} "
RS-9900A	Micro Sterilizing adjustable cushions 8" x 8 ^{1/2} " x 1 ^{1/4} "
RS-9902	Micro Sterilizing 8" x 12 ^{1/2} " x 1 ^{1/4} "
RS-9902A	Micro Sterilizing adjustable cushions 8" x 12 ^{1/2} " x 1 ^{1/4} "

An effective way of protecting delicate and micro surgical instruments for sterilization.

- Separating pins provide for instrument stability
- Silicone rubber inserts protect instruments
- Removable cover



RS-9900A Shown
Not Actual Size

RS-9910	Sentry Instrument Tray 9 ^{3/8} " x 8 ^{5/8} " x 1 ^{1/2} "
----------------	--

Protect your delicate instruments with the cost effective Sentry Instrument Tray. The Sentry Instrument Tray is a light-weight, inexpensive tray to autoclave all your instruments. The tray has a separate top and bottom that lock conveniently together. The tray is also lined with silicone extrusions to secure your instruments.

Stackable trays make it easy to use multiple trays at once within an autoclave.

The Sentry tray size is 9^{3/8}" x 8^{5/8}" x 1^{1/2}" and is made of Ultem 2000 Plastic.

