

# Clips and Clamps



		<b>RS-5420</b>	Micro Clip straight 70gr pressure clip width 0.75mm jaw length 5mm
		<b>RS-5422</b>	Micro Clip straight 70gr pressure clip width 0.75mm jaw length 6mm
		<b>RS-5424</b>	Micro Clip straight 70gr pressure clip width 1mm jaw length 8mm
		<b>RS-5426</b>	Micro Clip straight 70gr pressure clip width 1.5mm jaw length 10mm
		<b>RS-5431</b>	Micro Clip slight curve 70gr pressure clip width 1mm jaw length 6mm
		<b>RS-5433</b>	Micro Clip strong curve 70gr pressure clip width 1mm jaw length 6mm
		<b>RS-5435</b>	Micro Clip 45° angle 70gr pressure clip width 1mm jaw length 6mm


Please use clip applying forceps with these micro clips to avoid damage and for ease of use.



RS-5410  
Clip Applying Forceps  
for micro clips  
length 5<sup>1</sup>/<sub>2</sub>"




RS-5440  
Micro Clip  
applying forceps  
with lock  
length 5<sup>3</sup>/<sub>4</sub>"



RS-5450

SCHWARTZ clip applying forceps for temporary clips length 8"



75% of actual size

RS-5452

SCHWARTZ temporary clip straight, smooth 795gr pressure jaw width 1.7mm jaw length 8mm clip length 1"



RS-5457

SCHWARTZ temporary clip slight angle, smooth 795gr pressure jaw width 1.7mm jaw length 8mm clip length 1"



RS-5459

SCHWARTZ temporary clip strong angle, smooth 795gr pressure jaw width 1.7mm jaw length 8mm clip length 1"



**RS-6470**      Micro Vascular Clip  
 vessel diameter 1mm  
 occlusion pressure 30-35gr  
 length 8mm



**RS-6472**      Micro Vascular Clip  
 vessel diameter 1-2mm  
 occlusion pressure 35-40gr  
 length 11mm



**RS-6474**      Micro Vascular Clip  
 vessel diameter 2-3mm  
 occlusion pressure 50-60gr  
 length 14mm



**RS-6488**      Micro Vascular Approximator  
 with moveable clips  
 vessel diameter 1mm  
 occlusion pressure 30-35gr  
 length 8mm



**RS-6490**      Micro Vascular Approximator  
 with moveable clips  
 vessel diameter 1-2mm  
 occlusion pressure 35-40gr  
 length 11mm




**RS-6492**      Micro Vascular Approximator  
 with moveable clips  
 vessel diameter 2-3mm  
 occlusion pressure 50-60gr  
 length 14mm



**RS-6496**      Micro Clip Setting  
 Forceps  
 length 4<sup>3</sup>/<sub>4</sub>"



To avoid damage and for ease of use, please use the micro vascular clips with these clip setting forceps.





RS-5470

DIEFFENBACH  
 serrefine  
 straight, very delicate  
 300gr pressure  
 length 1 1/8"



RS-5471

DIEFFENBACH  
 serrefine  
 curved, very delicate  
 300gr pressure  
 length 1 1/8"




RS-7438

DIETRICH  
 micro bulldog clamp  
 straight  
 80gr pressure  
 jaw length 8mm  
 jaw width 1.2mm  
 length 2"



RS-7439

DIETRICH  
 micro bulldog clamp  
 curved  
 80gr pressure  
 jaw length 8mm  
 jaw width 1.2mm  
 length 2"



**RS-7440**      JOHNS HOPKINS  
 bulldog clamp  
 straight  
 700gr pressure  
 length 1<sup>1/2</sup>"



**RS-7441**      JOHNS HOPKINS  
 bulldog clamp  
 curved  
 700gr pressure  
 length 1<sup>1/2</sup>"



**RS-7442**      JOHNS HOPKINS  
 bulldog clamp  
 straight  
 800gr pressure  
 length 2"



**RS-7443**      JOHNS HOPKINS  
 bulldog clamp  
 curved  
 800gr pressure  
 length 2"



**RS-7444**      JOHNS HOPKINS  
 bulldog clamp  
 straight  
 800gr pressure  
 length 2<sup>1/2</sup>"




**RS-7445**      JOHNS HOPKINS  
 bulldog clamp  
 curved  
 800gr pressure  
 length 2<sup>1/2</sup>"



**RS-7446**      JOHNS HOPKINS  
 bulldog clamp  
 straight  
 1000gr pressure  
 length 3"



**RS-7447**      JOHNS HOPKINS  
 bulldog clamp  
 curved  
 1000gr pressure  
 length 3"





RS-7500

DeBAKEY  
atraumatic bulldog clamp  
straight  
400-500gr pressure  
jaw length 2cm  
overall length 7.5cm



RS-7511

DeBAKEY  
atraumatic bulldog clamp  
curved  
400-500gr pressure  
jaw length 2cm  
overall length 7.5cm






RS-7530  
 GLOVER  
 atraumatic bulldog clamp  
 straight  
 tension adjusting screw  
 240-260gr pressure  
 jaw length 2cm  
 overall length 7cm




RS-7541  
 GLOVER  
 atraumatic bulldog clamp  
 curved  
 with tension adjusting screw  
 240-260gr pressure  
 jaw length 2cm  
 overall length 7cm





RS-7760

JONES  
towel clamp  
cross action  
light weight  
length 2<sup>1</sup>/<sub>4</sub>"



RS-7784 Shown

RS-7780

RS-7784

BACKHAUS  
towel clamp  
length 3<sup>1</sup>/<sub>2</sub>"  
length 5<sup>1</sup>/<sub>4</sub>"

