



**GEMINI™**  
Cautery System

- RS-300** Gemini Two Temperature Cautery Kit contains: handle and 2 - 5mm tips and 2 - 10mm tips
- RS-301** 5mm replacement tips pack of 10
- RS-302** 10mm replacement tips pack of 10

**TECHNICAL DATA:**

Length (without tip): 160.75mm

Weight (with two AA batteries): 78g

Electrical: 3 volts, 1.5 volts

Tip Sizes: 5x 1mm, 10x 1mm



The GEMINI CAUTERY SYSTEM is the perfect tool for your cauterizing needs. The unique tip sizes combined with two different operating temperatures allow for maximum versatility. The convenient battery spacer converts the hand piece from high temperature to low temperature.



The Gemini Cautery Unit has the following features:

- Simple tip design for easy replacement
- Uses standard AA batteries
- Easy screw cap for convenient battery replacement
- Hard plastic case with die cut foam for convenient storage
- Battery spacer for optional temperature operation
- One year warranty

