


		RS-5420	Micro Clip straight 70gr pressure clip width 0.75mm jaw length 5mm
		RS-5422	Micro Clip straight 70gr pressure clip width 0.75mm jaw length 6mm
		RS-5424	Micro Clip straight 70gr pressure clip width 1mm jaw length 8mm
		RS-5426	Micro Clip straight 70gr pressure clip width 1.5mm jaw length 10mm
		RS-5431	Micro Clip slight curve 70gr pressure clip width 1mm jaw length 6mm
		RS-5433	Micro Clip strong curve 70gr pressure clip width 1mm jaw length 6mm
		RS-5435	Micro Clip 45° angle 70gr pressure clip width 1mm jaw length 6mm

Please use clip applying forceps with these micro clips to avoid damage and for ease of use.



RS-5410 Clip Applying Forceps
for micro clips
length 5¹/₂"



RS-5440 Micro Clip
applying forceps
with lock
length 5³/₄"

