



Not Actual Size

**RC-5903** Mini-Clipper with no. 0000 blade  
 UL Listed. Operates on 120 volts AC, weight 26 oz. Insulated 8 ft. cord, double insulated unit.

**RC-5903-RB** Replacement blade for RC-5903



Not Actual Size

**RC-5904** Cordless/Rechargeable Clipper with no. 0000 blade  
 UL Listed. Operates on battery or cord. Comes with recharging unit and 8 ft. cord. weight 1.2 lbs.

**RC-5904-RB** Replacement blade for RC-5904



RS-5905 Shown  
 Not Actual Size

**RC-5905** A-5 Clipper with detachable blade with No. 10 blade  
**RC-5906** with No. 40 blade  
 UL Listed. Operates on 120 volts AC weight 26 oz. This clipper is double insulated and meets OSHA specifications. Weight 26ozs.

**RC-5910** **A-5 replacement blades** medium-general hair length 1/16" size 10

**RC-5911** for clipping surgical area hair length 1/100" size 40

**RC-5912** for clipping surgical area hair length 1/125" size 50

**RC-5920** Spray Disinfectant, 14 oz.

**RC-5921** Spray Lubricant, 16 oz.



