



Disposable Feeding Needles plastic ball & hub
malleable cannula, shipped in boxes of 5

Cat. # Gauge & Length

FN-9921 22G - 1.5" - straight

FN-9918 18G - 2" - straight

REUSEABLE FEEDING NEEDLES
Stainless Steel

<i>Cat. #</i>	<i>Gauge & Length</i>	<i>Ball Diameter</i>	<i>Cat. #</i>	<i>Gauge & Length</i>	<i>Ball Diameter</i>
STRAIGHT			CURVED		
FN-7900	24G - 1"	1.25mm	FN-7910	20G - 1.5"	2.25mm
FN-7901	22G - 1"	1.25mm	FN-7911	20G - 3"	2.25mm
FN-7920	22G - 1.5"	1.25mm	FN-7913	18G - 2"	2.25mm
FN-7921	20G - 1"	2.25mm	FN-7914	18G - 3"	2.25mm
FN-7903	20G - 3"	2.25mm	FN-7916	16G - 3"	3mm
FN-7905	18G - 2"	2.25mm	FN-7917	14G - 3"	4mm
FN-7906	18G - 3"	2.25mm			
FN-7908	16G - 3"	3mm			
FN-7909	14G - 3"	4mm			
FN-7944	20G - 1.5"	2.25mm Malleable			

