



RS-7500

DeBAKEY
atraumatic bulldog clamp
straight
400-500gr pressure
jaw length 2cm
overall length 7.5cm



RS-7511

DeBAKEY
atraumatic bulldog clamp
curved
400-500gr pressure
jaw length 2cm
overall length 7.5cm



RS-7530

GLOVER
atraumatic bulldog clamp
straight
tension adjusting screw
240-260gr pressure
jaw length 2cm
overall length 7cm



RS-7541

GLOVER
atraumatic bulldog clamp
curved
with tension adjusting screw
240-260gr pressure
jaw length 2cm
overall length 7cm

