



RS-5470

DIEFFENBACH  
 serrefine  
 straight, very delicate  
 300gr pressure  
 length 1<sup>1/8</sup>"



RS-5471

DIEFFENBACH  
 serrefine  
 curved, very delicate  
 300gr pressure  
 length 1<sup>1/8</sup>"



RS-7438

DIETRICH  
 micro bulldog clamp  
 straight  
 80gr pressure  
 jaw length 8mm  
 jaw width 1.2mm  
 length 2"



RS-7439

DIETRICH  
 micro bulldog clamp  
 curved  
 80gr pressure  
 jaw length 8mm  
 jaw width 1.2mm  
 length 2"

