

Retractors





RS-6500 AGRICOLA
 3x3 prongs
 depth 4mm
 spread 1"
 length 1 1/2"



RS-6502 GOLDSTEIN
 3x3 prongs
 depth 5mm
 spread 3/4"
 length 1"



RS-6510 Shown

RS-6510 ALM
 4x4 prongs
 depth 3mm
 spread 2 1/4"
 length 2 3/4"
RS-6514 depth 13mm
 spread 3"
 length 4"



RS-6521 Shown


RS-6521 Cross Action Retractor
 with ratchet
 spread 1"
 length 4"
 angled
 4x4 blunt prongs
RS-6524 straight
 4x4 blunt prongs




RS-6540 BARRAQUER
 wire blades
 spread 12mm
 blade width 10mm
 length 1 3/4"




RS-6550 MURDOCK
 non-magnetic, chrome
 spread 38mm
 blade width 15mm
 length 2¹/₄"



RS-6560 MCKINNEY RETRACTOR
 spread 2"
 blade width 15mm




RS-6561 COOK
 non-magnetic
 spread 1.5"
 blade width 11mm
 length 2"



RS-6570 CASTROVIEJO
 non-magnetic
 spread 35mm
 blade length 15mm
 blade depth 5mm
 length 3³/₄"




RS-6585 DAVIS
 3x3 blunt prongs
 spread 2.5"
 length 4¹/₂"



RS-6610
RS-6612

SENN
double end
length 6^{1/2}"
3 sharp prongs
3 blunt prongs



RS-6610 Shown

RS-6614

RAGNELL
double end
3mm x 5mm
5mm x 15mm
length 6"

**RS-6628**

Rake retractor
rigid shaft
length 6"
2 prongs, blunt



RS-6622 Shown

1-800-424-2984
www.roboz.com



RS-6630 Rake retractor
rigid shaft
length 6"
3 prongs, blunt




RS-6624 Shown

RS-6670 DESMARRES
length 5 1/2"
blade width 11mm
RS-6672 blade width 13mm





RS-6672 Shown

RS-8600 WEITLANER
2 x 3 prongs
spread 1 1/2"
length 4"
sharp
RS-8602 blunt



RS-8600 Shown





RS-8612 Shown

RS-8612WEITLANER
3 x 4 prongs, blunt
spread 1^{3/4}"
length 5^{1/2}"**RS-8614**spread 3^{1/2}"
length 6^{1/2}"**RS-8616**spread 4^{3/4}"
length 8"**1-800-424-2984**
www.roboz.com