



Scalpels






- RS-9801-10** Sterile Scalpel Blades
- RS-9801-11** interchangeable
- RS-9801-12** nos. 10, 11, 12, 12D, 14,
- RS-9801-12D** 15, 15C
- RS-9801-14** for use with handles
- RS-9801-15** #3 & #7
- RS-9801-15C** box of 100



- RS-9861-18** Sterile Scalpel Blades
- RS-9861-19** interchangeable
- RS-9861-20** nos. 18, 19, 20, 21, 22,
- RS-9861-21** 23, 24, 25A, 36
- RS-9861-22** for use with handle #4
- RS-9861-23** box of 100
- RS-9861-24**
- RS-9861-25A**
- RS-9861-36**

Please specify blade number

1-800-424-2984
www.roboz.com



RS-9843 Scalpel Handle, #3
solid
length 4"




RS-9847 Scalpel Handle, #7
solid
length 6 1/4"



RS-9884 Scalpel Handle, #4
solid
length 4 3/4"




1-800-424-2984
www.roboz.com



37-1350

Handle
round, knurled
length 13cm
diameter .5cm



37-4310

Handle
round, knurled
length 10cm
diameter .5cm



37-4312

Handle
round, knurled
length 7.5cm
diameter .5cm



37-5910

Needle blade
3.0mm cutting edge
box/3



37-5920

Needle blade
5.0mm cutting edge
box/3



37-6100

Miniature blade
1.5mm cutting edge
box/24



37-7513 Miniature scalpel
disposable blade
assembly w/o handle
blue
15° 3mm depth
box/6




For use with handle #37-1350, 37-4310, or 37-4312

37-7515 Miniature scalpel
disposable blade
assembly w/o handle
green
15° 5mm depth
box/6



For use with handle #37-1350, 37-4310, or 37-4312

1-800-424-2984
www.roboz.com



**37-7514**

Miniature 1 piece
disposable scalpel
blue
15° 3mm depth
box/6

**37-7516**

Miniature 1 piece
disposable scalpel
green
15° 5mm depth
box/6

**37-7546**


Miniature 1 piece
disposable scalpel
tan
45° 3mm depth
box/6

1-800-424-2984
www.roboz.com




RS-6456

CASTROVIEJO
blade breaker
with push-pull lock
length 5^{1/4"}



1-800-424-2984
www.roboz.com



COMPARISON OF MILLIMETERS TO INCHES

MM	INCHES	MM	INCHES
0.33	.013	6	.236
0.5	.020	6.33	.249
0.67	.026	6.5	.256
1	.039	6.67	.263
1.33	.052	7	.275
1.5	.059	7.33	.288
1.67	.066	7.5	.295
2	.079	7.67	.302
2.33	.092	8	.315
2.5	.098	8.33	.328
2.67	.105	8.5	.335
3	.118	8.67	.341
3.33	.131	9	.354
3.5	.138	9.33	.367
3.67	.144	9.5	.374
4	.157	9.67	.380
4.33	.170	10	.394
4.5	.177	10.33	.407
4.67	.184	10.67	.420
5	.197	11	.433
5.33	.210	11.33	.446
5.5	.216	11.67	.459
5.67	.223	12	.472

COMPARISON OF INCHES TO MILLIMETERS

INCHES	MM	INCHES	MM
1/32	0.79	17/32	13.49
1/16	1.59	9/16	14.29
3/32	2.38	19/32	15.08
1/8	3.18	5/8	15.88
5/32	3.97	21/32	16.67
3/16	4.76	11/16	17.46
7/32	5.56	23/32	18.26
1/4	6.35	3/4	19.05
9/32	7.14	25/32	19.84
5/16	7.94	13/16	20.64
11/32	8.73	27/32	21.43
3/8	9.53	7/8	22.22
13/32	10.32	29/32	23.02
7/16	11.11	15/16	23.81
15/32	11.91	31/32	24.61
1/2	12.7	1"	25.4

COMPARISON OF INCHES TO CENTIMETERS

INCHES	CM
1	2.54
2	5.08
3	7.62
4	10.16
5	12.7
6	15.24
7	17.78
8	20.32
9	22.86
10	25.40
11	27.94
12	30.48
13	33.02
14	35.56
15	38.10

COMPARISON OF CENTIMETERS TO INCHES

CM	INCHES
1	.3937
2	.7874
3	1.1811
4	1.5748
5	1.9685
6	2.3622
7	2.7559
8	3.1496
9	3.5433
10	3.937
20	7.8740
30	11.811
40	15.748

