

with 2x Magnifications (2.5x with Accessory Lenses) coated lenses for minimum loss of light

RS-6680

Binocular magnifying loupe mounted on spectacle frame that can be fitted with the wearer's usual corrective lenses. with carrying case

FEATURES:

- Light weight
- 30cm focal distance
- 10cm viewing area
- Elimination of error of telescopic parallax
- Magnification 1.8x without distortion of position
- Comfortable 12" working distance
- Adjustable pupillary distance
- Individually adjustable lenses
- Comfortable peripheral vision
- Loupe flips up-out of field of vision




75% of Actual Size

RS-6682

Magnifying Loupe only

FEATURES:

- Light weight
- 30cm focal distance
- 10cm viewing area
- Elimination of error of telescopic parallax
- Magnification 1.8x without distortion of position
- Comfortable 12" working distance
- Adjustable pupillary distance
- Individually adjustable lenses
- Comfortable peripheral vision
- Loupe flips up-out of field of vision

**RS-6684**

Accessory Clip-on lenses
2.5x magnification
17cm focal distance
4.5cm viewing area

