



Disposable Feeding Needles plastic ball & hub malleable cannula, shipped in boxes of 5

Cat. # Gauge & Length
 FN-9921 20G - 1.5" - straight
 FN-9918 18G - 2" - straight

REUSEABLE FEEDING NEEDLES
Stainless Steel

REUSEABLE FEEDING NEEDLES
Size Recommendations

CAT. #	GAUGE & LENGTH	BALL DIAMETER
STRAIGHT		
FN-7900	24G - 1"	1.25mm
FN-7901	22G - 1"	1.25mm
FN-7920	22G - 1.5"	1.25mm
FN-7921	20G - 1"	2.25mm
FN-7903	20G - 3"	2.25mm
FN-7905	18G - 2"	2.25mm
FN-7906	18G - 3"	2.25mm
FN-7908	16G - 3"	3mm
FN-7909	14G - 3"	4mm
FN-7944	20G - 1.5"	2.25mm Malleable
CURVED		
FN-7910	20G - 1.5"	2.25mm
FN-7911	20G - 3"	2.25mm
FN-7913	18G - 2"	2.25mm
FN-7914	18G - 3"	2.25mm
FN-7916	16G - 3"	3mm
FN-7917	14G - 3"	4mm

LABORATORY ANIMAL		RECOMMENDED STANDARD SIZES			
SPECIES	WT. RANGE IN GRAMS	GAUGE	LENGTH	BALL DIAMETER	SHAPE
MICE	To 14g	24G	1"	1.25mm	straight
	15-20g	22G	1",1.5"	1.25mm	straight
	20-25g	20G	1",1.5",3"	2.25mm	straight,curved
	25-30g	18G	1",1.5",2"	2.25mm	straight,curved
	30-35g	18G	2",3"	2.25mm	straight,curved
RATS	50-75g	20G	1",1.5"	2.25mm	straight
	75-120g	18G	1",1.5"	2.25mm	straight,curved
	100-200g	18G	2",3"	2.25mm	curved
		16G	2"	3mm	straight,curved
	150-300g	16G	3",4"	3mm	curved
		14G	3"	4mm	curved
	200-350g	13G	3"	4mm	straight
HAMSTERS	60-200g	18G	2"	2.25mm	curved
GUINEA PIGS	250-300g	18G	1.5",2"	2.25mm	curved
	350-450g	16G	3",4"	3mm	curved
PIGS	400-600g	14G	3"	4mm	curved
		13G	3"	4mm	straight
	RABBITS	1-3kg	16G	3",4"	3mm
PRIMATES	To 4.5kg	12G	3.5"	6.25mm	curved

